

D-Check Bio-Acoustic Simulator

The D-Check Bio-Acoustic Simulator performs dependable daily calibration checks of the SMART TONE Automatic Audiometer, and can be used for OSHA, MSHA, or DOT hearing programs. Because the D-Check is powered by the SMART TONE audiometer, there is no need to worry about battery life. This advanced design simplifies the daily check and offers greater reliability and repeatability.



Features

- Quick and accurate daily calibration check
- Uses no batteries – power supplied by the audiometer
- Color-coded ear cups reduce confusion and errors
- Table-top or test booth wall-mount